

Release Notes

Viedoc 4.7

Available for users on
2014-08-06

Release Notes

Viedoc 4.7

New functionality



List of new functionality that has been added to this release.

- Queries are visible to everyone until approved or rejected.
- Queries now have a study sequence number and an item sequence number.
- Any field can be confirmed as missing by a user with data entry permission. Missing confirmations shall be reviewed by user with raise query permission (+1 task).
- Query audit trail: All changes to queries are audit trailed with user and date & time and shown in query popup.
- Any form or study event to be reviewed will yield one task for the reviewer for each of clinical review, data review and SDV
- Query export
- Data review export
- Signing console
 - Eye icon (has user seen the form)
 - Navigate to and from forms
- Data review console

Release Notes

Viedoc 4.7

Updated functionality and bug fixes



List of new functionality that has been updated in this release.

- Form functionality:
 - The form save button is now hidden in edit mode until a change has been made to the form
 - Improved performance and scalability
 - Bug fixes
- Query workflow:
 - Queries can be raised on function update fields
 - Queries can be raised on items with no label.
 - Form fields with queries are highlighted with an orange border
 - The form stays open if any queries are generated when saving the form
 - Validation queries: raised by the system when an expression in the form evaluates to false. Data can be confirmed correct or validation queries will be closed automatically if the data is corrected.
 - Validation queries are triggered on changes to any fields included in the expression, including cross-form fields
 - Validation queries are closed and re-raised (if expression evaluates to false) upon expression execution.
 - Validation can now be set on function update fields
 - Query resolution is saved even if the form is submitted with required fields missing.

Release Notes

Viedoc 4.7

Updated functionality and bug fixes



List of new functionality that has been updated in this release.

- Lock icon will be shown on subject card if status is 100% and all forms and events are locked
- Task count will be shown as 9+ if the number of tasks is greater than 9
- Numeric function update fields will format the returned function value according to the precision set by the designer
- Changes in function update fields will automatically add a change reason (Function execution).
- Functions will be executed for hidden fields that become visible.
- Default values will be populated for hidden fields that become visible.
- Form review status is hidden in edit mode.
- Review (clinical review, data review and SDV) and Lock form and study event date

The following limitations that have been reported in earlier versions of Viedoc are now corrected in this release:

- Validation will only re-evaluate when a change is made to the field where it is attached to
- Functions will only re-execute when a change is made to a field;
 - That is included in the function expression, AND
 - That is in the same form as the field where the function is attached

Release Notes

Viedoc 4.7

Known limitations



List of limitations that have been found with this and earlier releases of Viedoc, and will be changed in future releases.

- The export preview feature is heavy to load and might cause the application to feel unresponsive
- Only one (selected) language can be imported from each CDISC ODM file
- Only one (selected) design can be imported from each CDISC ODM file
- When the study workflow is saved, the language in the design settings are lost and have to be reset

Release Notes

Viedoc 4.7

Signoff



Signature Confirmation by
the product owner.

Uppsala, Sweden 2014-07-30

Henrik Blombergsson

PCG Solutions